

Tapones de silicona, con dos agujeros

DIN 12871.

Número de agujeros:

Color:

Material:

Resistencia a la temperatura:

2

natural

Silicona

-60 ... 180 °C

Ø fondo mm	Ø superior mm	Altura mm	Ø agujeros mm	Estándar	Resistencia máx. a la temperatura °C	Resistencia mín. a la temperatura °C
14,0	18,0	20,0	3,0	DIN 12871	180	-60
17,0	22,0	25,0	4,0	DIN 12871	180	-60
18,0	24,0	30,0	4,0	DIN 12871	180	-60
21,0	27,0	30,0	5,0	DIN 12871	180	-60
23,0	29,0	30,0	5,0	DIN 12871	180	-60
26,0	32,0	30,0	6,0	DIN 12871	180	-60
26,0	32,0	30,0	7,0	DIN 12871	180	-60
29,0	35,0	30,0	6,0	DIN 12871	180	-60
31,0	38,0	35,0	6,0	DIN 12871	180	-60
31,0	38,0	35,0	7,0	DIN 12871	180	-60
36,0	44,0	40,0	5,0	DIN 12871	180	-60
36,0	44,0	40,0	8,0	DIN 12871	180	-60
47,0	55,0	40,0	7,0	DIN 12871	180	-60
50,5	59,5	45,0	9,0	DIN 12871	180	-60
50,5	59,5	45,0	11,5	DIN 12871	180	-60
56,0	65,0	45,0	11,5	DIN 12871	180	-60
64,5	75,5	55,0	9,0	DIN 12871	180	-60