

Exchangeable probes for PDF data logger EBI 310 TX

Type	Min. temperature °C	Max. temperature °C	Temperature	Ø mm	Length mm	Dimensions (ØxL)	Cable length m	PK
TPX 310-P1	-200	200	-200 ... 200 °C	5	45	5 x 45 mm	3	1
TPX 310-P2	-50	180	-50 ... 180 °C	3	130	3 x 130 mm	3	1
TPX 310-P3	-50	180	-50 ... 180 °C	3	130	3 x 130 mm	1	1
TPX 310-P4	-100	400	-100 ... 400 °C	1.5	50	1.5 x 50 mm	3	1