

Viscometers ViscoQC 100

- 3.5" colour LC display
- Direct print / export

Type	Viscosity range mPas	Max. torque Ncm	PK	Cat. No.
ViscoQC 100 - L	15 ... 6000000	0.00673	1	6.313 808
ViscoQC 100 - R	100 ... 40000000	0.07187	1	6.313 809
ViscoQC 100 - H	800 ... 320000000	0.57496	1	6.313 810

Zubehör:

- Accessories for Viscometers ViscoQC / Anton Paar Germany
- Cone/plate measuring system PTD 100 Cone-Plate / Anton Paar Germany
- Temperature chamber ETD 300 for rotational viscometers VisQC / Anton Paar Germany

Alternativen:

- Viscometers ViscoQC 300 / Anton Paar Germany