

Pre-treatment cartridges Progard®; for water purification systems

The pre-treatment cartridges contain polyphosphates, activated carbon and depth filter for water pre-treatment and serve to protect reverse osmosis membrane against organic fouling, chlorine oxidation, and mineral scaling.

Type	For	PK	Cat. No.
Progard®; TS2	Elix®; Reference/Essential/Advantage, RiOs®; Essential, AFS®, Milli-Q®, Integral	1	6.253 376
Progard®; TL1-C	Elix®; 20/35/70/100, RiOs®; 30/50/100/150/200	1	6.224 448
Progard®; T3	Milli-Q®, Direct	1	7.970 086
Progard®; S2	Direct-Q®, Elix®, RiOs®, 5/8/16	1	6.240 434