

# Burettes, straight stopcock, brown glass

SILBERBRAND, Boro 3.3. Calibrated to deliver (TD, Ex). Schellbach stripe.

25 ml: Error limit according to Class B, DIN EN ISO 385.

50 ml: Class B, DIN EN ISO 385.

Accuracy class:

B  
Straight

Nominal capacity ml	Division ml	Tolerance ±·ml	Standard	Material stopcock	Colour	Material	With Schellbach stripe	PK	Cat. N
25	0.10	0.08	DIN EN ISO 385	PTFE key	Brown	Boro 3.3	Yes	2	4.689
50	0.10	0.10	DIN EN ISO 385	PTFE key	Brown	Boro 3.3	Yes	2	4.689