

LLG-Filter papers, qualitative, circles, slow

Pore size 2 to 4 µm.

Thickness:	0.16 mm
Filtration speed:	200
Particle retention:	2...4 µm
Ash content:	0.060 %

Ø mm	Weight g/m ²	Max. filtration speed s	Material	Type of analysis	Shape	PK	Cat. No.
70.0	80	200	Cellulose	Qualitative	Round filters	100	7.970 267
90.0	80	200	Cellulose	Qualitative	Round filters	100	7.970 268
100.0	80	200	Cellulose	Qualitative	Round filters	100	6.326 360
110.0	80	200	Cellulose	Qualitative	Round filters	100	6.242 672
125.0	80	200	Cellulose	Qualitative	Round filters	100	7.970 134
150.0	80	200	Cellulose	Qualitative	Round filters	100	7.970 139
185.0	80	200	Cellulose	Qualitative	Round filters	100	6.242 671
240.0	80	200	Cellulose	Qualitative	Round filters	100	7.971 615

Alternativen:

LLG-Filter papers, qualitative, circles, very fast / LLG Labware

LLG-Filter papers, qualitative, circles, medium fast / LLG Labware