

Embedding cassettes fine pore, pore style small holes, with attached lid

For very small biopsy specimens. The embedding cassettes are comprehensively validated for print quality, processing quality and reliability during embedding and cutting.

- Opening and closing forces have been carefully adjusted
- PrintMate AS printer compatible tubes are precisely engineered to reduce jamming and prevent downtime
- The colour spectrum was selected for reliability during scanning
- Lid hinge at the front edge

Dimensions (WxDxH): 28.5 x 41.2 x 6 mm

Material	Colour	Width mm	Depth mm	Height mm	Package	PK	Cat. No.
Plastic	White	28.5	41.2	6.0	1 x Bag of 1,000	1000	4.688 432
Plastic	Blue	28.5	41.2	6.0	1 x Bag of 1,000	1000	4.688 429
Plastic	Green	28.5	41.2	6.0	1 x Bag of 1,000	1000	4.688 430
Plastic	Pink	28.5	41.2	6.0	1 x Bag of 1,000	1000	4.688 431
Plastic	White	28.5	41.2	6.0	10 x Tube of 75	750	4.688 484
Plastic	Blue	28.5	41.2	6.0	10 x Tube of 75	750	4.688 481
Plastic	Green	28.5	41.2	6.0	10 x Tube of 75	750	4.688 482
Plastic	Pink	28.5	41.2	6.0	10 x Tube of 75	750	4.688 483