

Benchtop timer PRISMA 200

ABS plastic case, large, easy-to-read LCD display, stopwatch function, normal time display, selectable time resolution 999 mins., 99/100 mins. or 9 hrs., 59 mins., 59 secs., 1/10 sec., loud audible alarm, 2 channels with start/stop operation, either separately or together. Battery Type: AA (R6).

Dimensions (WxDxH): 130 x 95 x 175 mm
Colour: White
Weight: 390 g
Battery/rechargeable battery: 1 x 1.5 V, AA

Type	No. of channels	Display	Type of display	Width mm	Depth mm	Height mm	Type of auxiliary energy	PK	Cat. No.
PRISMA 200	2	Digital	LCD	130	95	175	Battery	1	9.260 178

Alternativen:

Stopclock, analogue, Mesotron / Hanhart 1882

Benchtop timer PRISMA 400 / Hanhart 1882