

# Embedding cassettes HistoScreen®; pore style mesh, with loose lid

For very small biopsy specimens. The embedding cassettes are comprehensively validated for print quality, processing quality and reliability during embedding and cutting.

- Opening and closing forces have been carefully adjusted
- PrintMate AS printer compatible tubes are precisely engineered to reduce jamming and prevent downtime
- The colour spectrum was selected for reliability during scanning
- Lid hinge at the back edge
- No sponges or bags necessary

Dimensions (WxDxH): 28.3 x 41 x 7 mm

Material	Colour	Width mm	Depth mm	Height mm	Package	PK	Cat. No.
Plastic	White	28.3	41.0	7.0	1 x Bag of 1,000	1000	4.688 412
Plastic	Blue	28.3	41.0	7.0	1 x Bag of 1,000	1000	4.688 409
Plastic	Green	28.3	41.0	7.0	1 x Bag of 1,000	1000	4.688 410
Plastic	Pink	28.3	41.0	7.0	1 x Bag of 1,000	1000	4.688 411
Plastic	Purple	28.3	41.0	7.0	1 x Bag of 1,000	1000	6.325 405
Plastic	White	28.3	41.0	7.0	10 x Tube of 75	750	4.688 452
Plastic	Blue	28.3	41.0	7.0	10 x Tube of 75	750	4.688 449
Plastic	Green	28.3	41.0	7.0	10 x Tube of 75	750	4.688 450
Plastic	Pink	28.3	41.0	7.0	10 x Tube of 75	750	4.688 451