

F96 MicroWell™; plates, PS, Nunclon™; Δ surface

Plates with excellent optical quality and flat well bottom. The Nunclon™; Δ cell culture-treated surface is hydrophilic and promotes cell adhesion and cell growth.

Material: PS
Bottom shape: Flat
Sterile: Yes

No. of wells	Nominal capacity µl	Colour	With lid	Package	PK	Cat. No.
96	400	Transparent	Yes	16 x Pack of 10	160	6.803 751
96	400	Transparent	Yes	Single packed	50	9.407 473
96	400	White	Yes	Single packed	50	6.802 646
96	400	White	Yes	16 x Pack of 10	160	6.238 692
96	400	Black	Yes	Single packed	50	6.802 647
96	400	Black	Yes	16 x Pack of 10	160	6.901 657
96	400	Transparent	No	18 x Pack of 10	180	4.008 767
96	400	Transparent	No	Single packed	50	9.407 475

Alternativen:

F96 MicroWell™ plates, PS, non-treated surface / Thermo Elect.LED

F96 MicroWell™ plates, PS, non-treated surface, with barcode / Thermo Elect.LED

F96 MicroWell™ plates, PS, treated surface / Thermo Elect.LED