

# Homogeniser set T 18 mini digital ULTRA-TURRAX<sup>&reg;</sup> Solution

**Scope of supply:** Homogeniser, dispersing head S 18 N - 19 G, stand R 1826,

boss head clamp R 182

Working volume:	1 ... 1,500 ml
Speed:	3,000 ... 25,000 1/min
Max. Viscosity:	2000 mPa·s
Noise level:	75 db(A)
Display:	Digital
With timer:	Yes
Dimensions (WxDxH):	66 x 100 x 186 mm
Weight:	1.8 kg
Rated capacity:	210 W
Power supply:	100 ... 240 V, 50/60 Hz
IP code:	IP 30

Type	Min. working volume ml	Max. working volume ml	Min. speed 1/min	M
T 18 mini digital ULTRA-TURRAX <sup>&amp;reg;</sup> Solution	1	1500	3000	25
T 18 mini digital ULTRA-TURRAX <sup>&amp;reg;</sup> Solution	1	1500	3000	25
T 18 mini digital ULTRA-TURRAX <sup>&amp;reg;</sup> Solution	1	1500	3000	25
T 18 mini digital ULTRA-TURRAX <sup>&amp;reg;</sup> Solution	1	1500	3000	25

## Zubehör:

Plate stands for overhead stirrers and Dispersers T 18 and T 25 / IKA-Werke  
Bossheads R 182 for overhead stirrers and Dispersers T 18 and T 25 / IKA-Werke  
Dispersing heads for T 18 ULTRA-TURRAX<sup>&reg;</sup> series / IKA-Werke  
Strap clamps RH 3 for overhead stirrers and Dispersers T 18 and T 25 / IKA-Werke

## Alternativen:

Homogeniser T 18 mini digital ULTRA-TURRAX<sup>&reg;</sup> / IKA-Werke