

Vacuum pump trolley, CP

Vacuum pump trolley with two condensate traps that can be operated alternately or together. The pump trolley consists of an aluminium frame mounted on lockable castors with a lower PE-shelf for a vacuum pump and a PE-top shelf which accommodates two cold traps. The condensate traps can each hold approx. 150ml. Spherical ground joints are sealed with FEP covered O-rings. Supplied without pump.

The standard assembly is designed for use with LN₂ as coolant. When using CO₂ and acetone please order 2 x CO₂ grid (Cat. No. 9.032 081).

Note: Price does not include vacuum pump.

Type	Description	PK	Cat. No.
CP 1	With pressure gauge	1	9.881 380
CP 2	Without pressure gauge	1	9.881 381