

Immersion thermostat Universa PRO

Temperature: 30 ... 200 °C
Temperature stability: ±0,02 K
Max. flow rate: 22 l/min
Max. delivery pressure: 0,550 bar
Type of display: TFT
Dimensions (WxDxH): 195 x 234 x 333 mm
Weight: 4.1 kg
Heating capacity: 2.8 kW
Plug type: CEE
Power supply: 200 ... 240 V, 50/60 Hz
IP code: IP 21

| Type | Min. temperature °C | Max. temperature °C | Programmable | Width mm | Depth mm | Height mm | Type of auxiliary energy connection | P |
|--------------|------------------------|------------------------|--------------|-------------|-------------|--------------|---|---|
| Universa PRO | 30 | 200 | Yes | 195 | 234 | 333 | Mains connection | 1 |

Alternativen:
Immersion thermostat Universa PRO, WLAN capable / Lauda Dr.R.Wobser