

# Condensers acc. to Dimroth for extractors, with PP olives, DURAN<sup>&reg;</sup> tubing

Made of DURAN<sup>&reg;</sup> tubing. For extractors and for condensing vapours.  
With cone. Water connection with thread GL14, with two threaded PP olives.  
Connection: PP olive

Ground size	Cooling surface cm <sup>2</sup>	Material	Type of construction	PK	Cat. No.
NS29/32	140	DURAN <sup>&amp;reg;</sup> tubing	Dimroth condenser	1	9.043 021
NS34/35	160	DURAN <sup>&amp;reg;</sup> tubing	Dimroth condenser	1	9.043 022
NS45/40	370	DURAN <sup>&amp;reg;</sup> tubing	Dimroth condenser	1	9.043 023
NS60/46	680	DURAN <sup>&amp;reg;</sup> tubing	Dimroth condenser	1	9.043 026
NS71/51	930	DURAN <sup>&amp;reg;</sup> tubing	Dimroth condenser	1	9.043 027

## Zubehör:

PTFE sealing rings for jointed glassware / Dr. Dietmar Glindemann  
Sleeves, cone joint, PTFE / Bohlender  
Sleeves with ribs, PTFE / Bohlender  
Joint clips KECK, POM, for conical ground joints / DWK Life Sciences  
Joint clips KECK, POM, for spherical ground joints / DWK Life Sciences  
Ground joint clips for sleeve connections, wire, Chrome-nickel steel / H W S Hans W. Schmidt  
Joint clips for conical joint sleeve connections, stainless steel, with set screw / Lenz-Laborglas  
Joint grease, K.W.S / J.P. PÅllath  
BUDDE - joint grease, silicone-free / Buddeberg  
Engler distilling flasks, DURAN<sup>&reg;</sup>, ASTM, with beaded rim / DWK Life Sciences  
Gas leak detection fluid, bottle / Otto Fahsig

## Alternativen:

Dimroth condensers for extractors, with glass olives, DURAN<sup>&reg;</sup> tubing / Lenz-Laborglas  
Condensers acc. to Allihn with GL threads, DURAN<sup>&reg;</sup> tubing / Lenz-Laborglas  
Condensers acc. to Liebig, DURAN<sup>&reg;</sup> tubing, with plastic olives / Lenz-Laborglas