

# Multi channel pipette mLine<sup>&reg;</sup>, mechanical, variabel

The mechanical pipettes of Sartorius are known for their ergonomic, user-friendly design and reliability. mLine<sup>&reg;</sup> pipettes represent an ergonomic choice for quality-conscious users who perform long lasting and repetitive manual pipetting. These pipettes are used in (bio)pharmaceutical laboratories, research institutes, universities, health care and industrial laboratories for various liquid handling tasks.

Capacity	Min. capacity $\mu\text{l}$	Max. capacity $\mu\text{l}$	Division $\mu\text{l}$	Accuracy $\pm\text{-R}\%$	No. of channels	PK	Cat. No.
1 ... 10 $\mu\text{l}$	1.0	10.0	0.010	1,5	8	1	6.262 705
1 ... 10 $\mu\text{l}$	1.0	10.0	0.010	1,5	12	1	6.241 712
10 ... 100 $\mu\text{l}$	10.0	100.0	0.100	0,9	8	1	7.660 044
30 ... 300 $\mu\text{l}$	30.0	300.0	0.200	0,6	8	1	6.228 996
30 ... 300 $\mu\text{l}$	30.0	300.0	0.200	0,6	12	1	6.234 097

Zubehör:  
Holder for Microliter pipettes / Sartorius Lab Instruments(Bio)