

Refill kits Salvequick® plaster strips

SALVEQUICK® plaster strips, waterproof - Refill 6036

When pulling out, one half of the plaster is ready to be used and must be consumed immediately. 45 adhesive bandages per pack: 27 strips of 72 x 19 mm, 18 strips of 72 x 25 mm

SALVEQUICK® plaster strips, elastic - Refill 6444

When pulling out, one half of the plaster is ready to be used and must be consumed immediately. 40 adhesive bandages per pack: 24 strips of 72 x 19 mm, 16 strips of 72 x 25 mm

SALVEQUICK® patch strips XL, elastic - Refill 6470

When pulling out, one half of the plaster is ready to be used and must be consumed immediately. 21 adhesive bandages per pack: 14 strips of 80 x 30 mm, 7 strips of 80 x 60 mm

SALVEQUICK® sterile plaster strips, detectable - Refill 6735CAP

When pulling out, one half of the plaster is ready to be used and must be consumed immediately. 35 sterile adhesive bandages per pack: 21 strips 72 x 19 mm, 14 strips 72 x 25 mm

Description	For	Sterile	Detectable	PK	Cat. No.
plaster strips, waterproof	SALVEQUICK®	No	No	45	6.902 277
plaster strips, elastic	SALVEQUICK®	No	No	40	7.658 508
patch strips, elastic	SALVEQUICK®	No	No	21	7.629 353
sterile plaster strips, detectable	SALVEQUICK®	Yes	Yes	35	7.622 837