

Smart vertical autoclaves TLV Basic

Temperature: 100 ... 140 °C
Max. pressure: 2.6 bar
With vacuum facility: Yes
With recooling: No
Dim. Int.(Ø x H): 400 x 650 mm
Ext. dim. (WxDxH): 610 x 870 x 1,110 mm
Weight: 151 kg
Heating capacity: 3600 W
Power supply: 230 V, 50/60 Hz

Type	Nominal capacity l	Min. temperature °C	Max. temperature °C	Int. Ø mm	Depth int. mm	Ext. width mm	Ext. depth mm	Ext. height mm
TLV-75	83	100	140	400	650	610	870	1110

Zubehör:
Thermoprinter IT/TLV, for automatic autoclaves / RT Espinal S.L. (Raypa)